

**Efficient, safe, and cost-effective
production of potent AAV vectors
using whole plants as bioreactors.**

Daniel Gibbs Ph.D.



Gene Therapy's Manufacturing Demand Challenge



If (GTx manufacturers) can't **streamline manufacturing**.....gene therapies will become **increasingly inaccessible** to the **patients that rely on them.**"

- FDA CBER Director Peter Marks

Example - Cost of 1 dose of Lenmeldy - \$4.25 million and Hemgenix - \$3.5 million
Cost of manufacturing 1 dose of Hemgenix is ~1-2 million

DeGroot, L. (2023, November 6). "Pull every lever": Marks doubles down on urgency to improve gene therapy manufacturing. Endpoints News. <https://endpts.com/pull-every-levermarks-doubles-down-on-urgency-to-improve-gene-therapy-manufacturing/>

Hosseini, M., Fath, S. ., Hosseini, M., & Fath, S. . (2023, July 6). Cutting the cost of gene therapy manufacturing. Roland Berger. <https://www.rolandberger.com/en/Insights/Publications/Cutting-the-cost-of-genetherapy-manufacturing.html>



Gene Therapy's Manufacturing Demand Challenge

2023 Bioreactor Capacity Demand Soared to 2 Billion Liters,
While the Average Batch Size Remained at Just 500L



2 Billion Liters = ~800 Olympic
swimming pools




500L
Bioreactor



Gene Therapy's Manufacturing Demand Challenge

To produce 200L Bioreactor run. $\sim 2e16$ vg total yield ($\sim 1e14$ vg/L)


New AAV Process Development



Lyle et al., 2023 *Biotechnology Bioengineering* 2023;1-14.
DOI: 10.1002/bit.28402



Cirsium's Plant-Based Solutions




With cutting-edge production technology, we utilize transient gene delivery in plants and high-density Controlled Environment Agriculture (CEA) to manufacture high-quality AAV vectors for genetic medicine.



Cirsium's Plant-Based Solutions

Produce $\sim 2e16$ vg of GMP Compatible Crude – eq. 200L Bioreactor Run ($\sim 1e14$ vg/L)




Plant-Produced AAVs Have Comparable Quality to HEK293T Cell-Produced AAVs

Properties	Plant-AAV9-CMV-GFP	Mammalian-AAV9-CMV-GFP	Validation
% Full particles in the crude lysate	31.87%	29.62%	ELISA (in house)
Genomic titer (vg/mL) (post purification)	2.5e12 vg/mL (qPCR) 2.5e12 vg/mL (dPCR)	2.5e12 vg/mL (qPCR) 2.6e12 vg/mL (dPCR)	qPCR (in house) dPCR(outsourced)
Endotoxin levels <USP 85>	<0.05EU/mL (Pass)	<0.05EU/mL (Pass)	LAL chromogenic test (outsourced)
Bacteriostasis and Fungistasis <USP 71>	Successfully validated for Sterility testing	Successfully validated for Sterility testing	Direct Transfer (outsourced)
Sterility Test <USP 71>	No growth observed (Pass)	No growth observed (Pass)	Direct Transfer (outsourced)
Mycoplasma Test	No mycoplasma detected	No mycoplasma detected	PCR based (in house) Culture (outsourced)




High titer/quality AAV vectors can be produced by using plant based AAV manufacturing


HEK293-Produced
AAV2-CMV-EGFP



Plant-Produced
AAV2-CMV-EGFP




↗ = full ↘ = empty



TEM Processing and
Imaging performed by
Salk Biophotonics
Core Facility


Plant-AAV9



Results

Aggregate intensity	3.4	%
Capsid intensity	96.6	%
Small particle intensity	0	%
Capsid diameter	30.4	nm
Total capsid titer	2.56E+12	cp/ml

HEK-AAV9



Results

Aggregate intensity	0.2	%
Capsid intensity	99.8	%
Small particle intensity	0	%
Capsid diameter	27.0	nm
Total capsid titer	2.64E+12	cp/ml

Unchained Labs Stunner data:

- PDI is <0.1 for both material, suggesting no aggregation
- DLS shows comparable capsid diameter
- Capsid titer - ~2.5e12/mL

Study design for comparability testing of plant-AAV with HEK293-AAV in non-human primates (NHP)

Table 1: Animal Recruitment Summary

Groups: Vehicle (n=1)
Plant-AAV (n=3)
HEK293-AAV (n=3)

RoA: Intravenous (2.5e12vg/kg)
Intrathecal (2.5e12vg)


In-life: 45days

All animals had no/low AAV9-Nab at the time of dosing.

Group	Animal ID	Treatment	Route	Dose Volume	Dose (vg/kg)	Baseline Body Weight (kg)
1	E128	Vehicle	IT	1 mL/kg	2.5 x 10 ¹²	5.61
			IV	1 mL		
2	E129	PLANT-AAV9-CMV-EGFP	IT	1 mL/kg	2.5 x 10 ¹²	5.44
			IV	1 mL		
	E153		IT	1 mL/kg	2.5 x 10 ¹²	6.34
			IV	1 mL		
	E352		IT	1 mL/kg	2.5 x 10 ¹²	6.11
			IV	1 mL		
3	E402	HEK293-AAV9-CMV-EGFP	IT	1 mL/kg	2.5 x 10 ¹²	5.33
			IV	1 mL		
	E409		IT	1 mL/kg	2.5 x 10 ¹²	5.63
			IV	1 mL		
	E413		IT	1 mL/kg	2.5 x 10 ¹²	6.04
			IV	1 mL		



AAV9-CMV GFP mRNA expression data in NHP tissue by qRT-PCR: fold change of plant AAV9 to HEK AAV9



- Quantification of GFP expression across tissue samples indicates higher transgene expression for plant AAV9 vs HEK AAV9 in most skeletal muscles and lumbar spinal cord, lymph nodes and testes.





Increased/comparable production of GFP Protein by plant-AAV relative to HEK293-AAV in the lumbar spinal cord tissue (NHP)



Tissue: Lumbar Spinal Cord

Groups: Vehicle (n=1)
Plant-AAV (n=3)
HEK293-AAV (n=3)

In-life: 45days



Assessment of inflammatory markers GFAP and ALT in NHP tissues demonstrates excellent tolerability of plant AAV9




Quantification of relative GFAP expression, a marker of glial reactivity, indicates no significant difference in both treatment groups from vehicle only control group

Quantification of relative ALT expression, a marker of inflammatory liver damage, also indicated no significant difference in both treatment groups compared to vehicle only

Conclusion: Both plant and HEK produced AAV9-CMV-GFP vectors were well tolerated in the CNS and liver in treated NHPs



Neutralizing Ab data (Luciferase assay) indicates significantly reduced antibody response to Plant-AAVs, compared to the HEK-AAVs



Summary Slide: Plants are the Answer to Gene Manufacturing

- Manufacturing AAVs in plants can reduce the cost of manufacturing by >10X.
- Since plant produced AAVs are linearly scalable, it can significantly reduce the process development times.
- Since downstream processing for plant produced AAVs is identical to existing mammalian DSP, this platform can be easily adopted.
- Plant produced AAVs are high quality and pass the QC requirements.
- Plant produced AAVs can transduce the correct cell types in non-human primates and show comparable biodistribution as that of mammalian cell produced AAVs.




Contact

Daniel Gibbs Ph.D.

Founder and CEO

dgibbs@cirsiumbio.com

**Recipient of ARPA-H BAA Award:
VECTORS Project. 140D042490002**



Visit

<https://Cirsiumbio.com>

